



TRANSFORMTM – A CUSTOM SHAPE CREATED BY YOU

Heat Mouldable Impression Trays



- Save time & money
- Helps overcome the clinical challenge of taking impressions where unusual soft tissue, bone shape or tooth position exists
- Mould chairside in less than 60 seconds
- Single Use
- Comfortable for patient
- Simple to use – dip in hot water at 70°C for 20 seconds and the material becomes soft and malleable
- Central pad allows use of finger or instrument to apply pressure
- Good retention slots
- Integral handle with thumb pad

Available in 3 colour-coded sizes in dentate and edentulous shapes



Item #	Description	
118201	Transform Tray - Dentate - Large Lower #1	12 pcs.
118202	Transform Tray - Dentate - Large Upper #2	12 pcs.
118203	Transform Tray - Dentate - Medium Lower #3	12 pcs.
118204	Transform Tray - Dentate - Medium Upper #4	12 pcs.
118205	Transform Tray - Dentate - Small Lower #5	12 pcs.
118206	Transform Tray - Dentate - Small Upper #6	12 pcs.
118207	Transform Tray - Edentulous - Large Lower #7	12 pcs.

Item #	Description	
118208	Transform Tray - Edentulous - Large Upper #8	12 pcs.
118209	Transform Tray - Edentulous - Medium Lower #9	12 pcs.
118210	Transform Tray - Edentulous - Medium Upper #10	12 pcs.
118211	Transform Tray - Edentulous - Small Lower #11	12 pcs.
118212	Transform Tray - Edentulous - Small Upper #12	12 pcs.
118213	Transform Tray - Edentulous - Mixed Bag 2 x #7-12	12 pcs.
118214	Transform Tray - Dentate - Mixed Bag 2 x #1-6	12 pcs.